

Code No.22-0002

## SPECIFICATION

## Isovaleraldehyde

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Density (20°C)	0.795~0.802g/ml
Refractive index 20°C	1.385~1.391
Water	max.0.5%
Acidity [as (CH <sub>3</sub> ) <sub>2</sub> CHCH <sub>2</sub> COOH]	max.1.5%
Assay (capillary GC)	min.95.0%