

Code No.22-0094

SPECIFICATION

Vinyl Cinnamate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ pale yellow, clear ~ nearly clear liquid to pass test
Solubility in ethanol	
Density (20°C)	1.067~1.074g/ml
Refractive index 20°C	1.583~1.589
Water	max.0.2%
Acidity (as C ₆ H ₅ CH:CHCOOH)	max.1.0%
Assay (capillary GC)	min.98.0%