

## SPECIFICATION

## Wakosil(R) 75

for Column Chromatography

REQUIREMENT	SPECIFICATION
Appearance	White, small granules or powder
Particle size (45~105 $\mu\text{m}$ )	min.70%
Particle size (coarser than 105 $\mu\text{m}$ )	max.20%
Particle size (finer than 45 $\mu\text{m}$ )	max.10%
Precipitation volume	under investigation
Loss on drying at 180°C	under investigation
pH (50g/l aqueous slurry,25°C)	under investigation
Flowing rate	under investigation
Separation power	to pass test
Activity	to pass test