

## SPECIFICATION

WakogelR 60N, 38~100  $\mu$ m

for Column Chromatography

REQUIREMENT	SPECIFICATION
Appearance	White powder
Particle size (38~100 $\mu$ m)	min.80%
Particle size (coarser than 100 $\mu$ m)	max.10%
Precipitation volume	2.2~3.0mL/g
Solubility in ethyl acetate	to pass test
Loss on drying at 180°C	max.7.0%
pH(50g/L aqueous slurry, 25°C)	6.0~7.5
Flowing rate	0.8~1.5mL/min
Separation power	to pass test
Activity	to pass test