

Code No.23-0277

## SPECIFICATION

WakosilR 60, 38~63  $\mu\text{m}$   
for Column Chromatography

REQUIREMENT	SPECIFICATION
Appearance	White small granules
Particle diameter (d50)	45~75 $\mu\text{m}$
Particle diameter (d10)	min.40 $\mu\text{m}$
Particle diameter (d90)	max.75 $\mu\text{m}$
Precipitation volume	1.9~2.5mL/g
Loss on drying at 105°C	4.0~10.0%
pH(50g/L aqueous slurry, 25°C)	6.5~7.5
Flowing rate	1.2~1.8mL/min
Separation power	to pass test
Activity	to pass test