

Code No.26-0167

SPECIFICATION

Zinc Stearate, 1.5 μ m

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly yellow, powder
Particle diameter (d50)	max.1.5 μ m
Loss on drying at 105°C	max.0.5%
Acidity [as CH ₃ (CH ₂) ₁₆ COOH]	max.0.5%
Assay (as Zn) (volumetric analysis)	10.0~11.0%
Assay (capillary GC)	under investigation