

| | | | | | | | |
|------------------------|------------------|---------|----------------|----------------------|-----------|----------------|------------|
| Frequency (MHz) | 399.76 | Nucleus | ¹ H | Number of Transients | 32 | Origin | JNM-ECS400 |
| Original Points Count | 65536 | Owner | delta | Points Count | 65536 | Pulse Sequence | proton.ho |
| Solvent | TETRACHLOROETHAN | | | Spectrum Offset (Hz) | 1998.9109 | Scan Rate (Hz) | 15024.94 |
| Temperature (degree C) | 30.000 | | | | | | |

