



DFIL0 Q1H  
 COMNT Biterbandol  
 EXMOD SGNON  
 OBNUC 1H  
 : OBFIN 5400.0 Hz  
 POINT 32768  
 FREQU 8103.7 Hz  
 SCANS 16  
 : ACQTM 2.022 sec  
 : PD 2.978 sec  
 : PM1 4.3 us  
 IRFIN 5400.0 Hz  
 IRATN 15  
 IRRPW 60 us  
 TEMP. 27.0 c  
 SLVNT CDCL3  
 EXREF 0.00 ppm  
 BF 0.25 Hz  
 RGAIN 17  
 XE 2701.6630 Hz  
 XS -227.2593 Hz



品名  
 ビテルタノール