


OBNUC
EXMOD
OBFRQ
OBSET
OBFIN
POINT
FREQU
SCANS
ACQTM
PD
PW1
IRNUC
CTEMP
SLVNT
EXREF
BF
RGAIN

1H
NON
399.65 MHz
124.00 KHz
10500.00 Hz
16384
7992.01 Hz
8
2.0500 sec
4.9500 sec
6.50 usec
1H
26.1 c
CDCL3
0.00 ppm
0.12 Hz
19



品名
クロフロップ