


1H  
proton.jpg  
399.78 MHz  
4.19 KHz  
7.29 Hz  
65535  
15024.04 Hz  
16  
4.3621 sec  
5.0000 sec  
4.90 usec  
1H  
24.1 c  
CDCL3  
0.00 ppm  
0.10 Hz  
50

OBNUC  
EXMOD  
OBFRQ  
OBSET  
OBFIN  
POINT  
FREQU  
SCANS  
ACQTM  
PD  
PWI  
IRNUC  
OTEMP  
SLVNT  
EXREF  
BF  
RGAIN



品名
フィブロニル