


1H  
 399.65 MHz  
 124.00 KHz  
 10500.0 Hz  
 32768  
 8000.0 Hz  
 8  
 4.096 sec  
 2.904 sec  
 6.2 us  
 25.3 c  
 DMSO  
 0.00 ppm  
 0.12 Hz  
 12

OBNUC 1H  
 EXMOD NON  
 OBFRQ  
 OBSET  
 OBF IN  
 POINT  
 FREQU  
 SCANS  
 ACQTM  
 PD  
 PW1  
 IRNUC 1H  
 CTEMP  
 SLVNT DMSO  
 EXREF  
 BF  
 RGAIN



品名
フェンブプロピモルフ