



OBNUC 1H  
EXMOD NON  
OBFRQ 399.65 MHz  
OBSET 124.00 KHz  
OBFIN 10500.0 Hz  
POINT 32768  
FREQU 8000.0 Hz  
SCANS 8  
ACQTM 4.096 sec  
PD 2.904 sec  
PW1 6.2 us  
IRNUC 1H  
CTEMP 26.5 c  
SLVNT DMSO  
EXREF 0.00 ppm  
BF 0.12 Hz  
RGAIN 13

品 名
(Z)-フェンプロキシメート