

|                         |                 |         |                |                      |           |                  |            |
|-------------------------|-----------------|---------|----------------|----------------------|-----------|------------------|------------|
| Frequency (MHz)         | 399.78          | Nucleus | <sup>1</sup> H | Number of Transients | 16        | Origin           | JNM-EC5400 |
| Original F1 File Count  | 65536           | Owner   | delta          | Points Count         | 65536     | File Sequence    | prcton.bp  |
| Solvent                 | METHYLENE-CHLOR |         |                | Spectrum Offset (Hz) | 1998.9109 | Since Shift (Hz) | 19024.04   |
| Temperature (degrees C) | 25.000          |         |                |                      |           |                  |            |

品 名  
インドキサカルブ-MP

