


ORNUC 1H  
 QBFIN 5400.0 Hz  
 P1INT 32768  
 FREQ 600.4 Hz  
 SCANS 8  
 ACQTM 2.730 sec  
 PD 5.000 sec  
 PW1 4.7 us  
 IRFIN 5400.0 Hz  
 INATN 1  
 IRRPW 50 us  
 TEMP. 27.0 C  
 SLVNT CDCL3  
 EXREF 0.00 ppm  
 BF 0.10 Hz  
 RGAIN 18  
 XE 2700.0000 Hz  
 XS -200.0000 Hz



品名
ナプロパミド