


OBNUC 1H  
EXMOD NON  
OBFRQ 399.65 MHz  
OBSET 124.00 KHz  
OBFIN 10500.00 Hz  
POINT 32768  
FREQU 7992.01 Hz  
SCANS 8  
ACQTM 4.1001 sec  
PD 4.9500 sec  
PW1 6.70 usec  
IRNUC 1H  
CTEMP 23.9 c  
SLVNT DMSO  
EXREF 0.00 ppm  
BF 0.12 Hz  
RGAIN 15



品名
プロペナゾール