


OBNUC 1H  
 OBFIN 5400.0 Hz  
 POINT 32768  
 FREQU 6002.4 Hz  
 SCANS 8  
 ACQTM 2.730 sec  
 PD 5.000 sec  
 PM1 6.0 us  
 IRFIN 5400.0 Hz  
 IRATN 0  
 IRRPW 50 us  
 TEMP. 27.0 c  
 SLVNT CDCL3  
 EXREF 0.00 ppm  
 BF 0.20 Hz  
 RGAIN 16  
 XE 2700.0000 Hz  
 XS -200.0000 Hz



品名
チオジカルブ