


1H
 proton.jpg
 399.78 MHz
 4.19 KHz
 7.29 Hz
 85535
 15024.04 Hz
 18
 4.3621 sec
 5.0000 sec
 4.90 usec
 1H
 24.5 c
 DMSO
 0.00 ppm
 0.12 Hz
 44

OBNUC
 EXMOD
 OBFRQ
 OBSET
 OBFTN
 POINT
 FREQU
 SCANS
 ACQTM
 PD
 PW1
 IRNUC
 CTEMP
 SLVNT
 EXREF
 BF
 RGAIN



品名
チアベンダゾール