


OBNUG 1H  
 EXMOD NON  
 OBFRQ 399.65 MHz  
 OBSET 124.00 KHz  
 OBF IN 10500.0 Hz  
 POINT 32768  
 FREQU 8000.0 Hz  
 SCANS 8  
 ACQTM 4.096 sec  
 PD 2.904 sec  
 PW1 6.2 us  
 IRNUC 1H  
 CTEMP 26.6 c  
 SLVNT CDCL3  
 EXREF 0.00 ppm  
 BF 0.12 Hz  
 RGAIN 15



品名
DMTPオキソン