


1H
 399.65 MHz
 124.00 KHz
 10500.0 Hz
 32768
 8000.0 Hz
 8
 4.096 sec
 2.904 sec
 6.2 us
 27.5 c
 7.00 ppm
 0.12 Hz
 14

OBRNUC 1H
 EXMOD NON
 OBFREQ
 OBSETE
 OBFIN
 POINT
 FREQU
 SCANS
 ACQTM
 PD
 PW1
 IRNUC 1H
 CTEMP
 SLVNT CDCL3
 EXREF
 BF
 RGAIN



品名
ジアフェンチウロン尿素