


1H
 ORNUC 1H
 EXMOD proton_jrp
 ORFREQ 399.78 MHz
 ORSET 4.19 KHz
 ORFIN 7.29 Hz
 POINT 16384
 FREQU 7503.00 Hz
 SCANS 8
 ACQTM 2.1837 sec
 PD 5.0000 sec
 PW1 5.00 msec
 1H 24.6 c
 DMSO
 SLVNT 6.00 ppm
 EXREF 0.12 Hz
 BF 40
 RGAIN



品名
トリチコナゾール