


1H proton.jpg
 399.78 MHz
 4.19 KHz
 7.29 Hz
 65535
 15024.04 Hz
 16
 4.3621 sec
 5.0000 sec
 4.90 usec
 1H 23.1 °
 DMSO 0.00 ppm
 0.12 Hz
 40

OBNUC
 EXMOD
 OBIFRG
 OBSET
 OBFIN
 POINT
 FREQU
 SCANS
 ACQTM
 PD
 PW1
 IRNUC
 CTEMP
 SLVNT
 EXREF
 BF
 RGAIN



品名
フェニピラザミン