


1H
 proton.jxp
 399.78 MHz
 4.19 KHz
 7.29 Hz
 65535
 15024.04 Hz
 16
 4.3621 sec
 5.0000 sec
 4.80 usec
 1H
 23.2 c
 DMSO
 2.50 ppm
 0.10 Hz
 42

OBNUC
 EXMOD
 OBFRQ
 OBSET
 OBFIN
 POINT
 FREQU
 SCANS
 ACQTM
 PD
 PW1
 IRNUC
 CTEMP
 SLVNT
 EXREF
 BF
 RGAIN



品名
フルオロピコリド