

| | | | | | | | |
|-----------------------|------------|----------------------|----------------|----------------------|----------|------------------------|------------------|
| Frequency (MHz) | 399.78 | Nucleus | ¹ H | Number of Transients | 8 | Origin | JNM-ECS400 |
| Original Points Count | 39979 | Owner | delta | Points Count | 65536 | Pulse Sequence | qNMR_experiment. |
| Solvent | ACETONE-d6 | Spectrum Offset (Hz) | 1998.9109 | Sweep Width (Hz) | 10017.28 | Temperature (degree C) | 25.000 |

| |
|------------|
| 品 名 |
| ダイファシノン標準品 |

