

| | | | | | | | |
|-----------------------|------------|----------------------|----------------|----------------------|----------|------------------------|------------|
| Frequency (MHz) | 399.78 | Nucleus | ¹ H | Number of Transients | 18 | Origin | JNM-EGS400 |
| Original Points Count | 65536 | Owner | delta | Points Count | 65536 | Pulse Sequence | proton.jxp |
| Solvent | ACETONE-d6 | Spectrum Offset (Hz) | 2010.1742 | Sweep Width (Hz) | 15024.04 | Temperature (degree C) | 25.000 |

品名
2,4-PA-イソプロピル

