


File Name

Frequency (MHz)	600.17	Nucleus	<sup>1</sup> H	Number of Transients	8	Origin	JNM-ECA600
Original Points Count	60018	Owner	delta	Points Count	65536	Pulse Sequence	qNMR_experiment.
Solvent	ACETONITRILE-d3	Spectrum Offset (Hz)	3014.5381	Sweep Width (Hz)	15007.50	Temperature (degree C)	25.000

品名
イプフルフェノキン



Frequency (MHz)	600.17	Nucleus	<sup>1</sup> H	Number of Transients	8	Origin	JNM-ECA600
Original Points Count	60018	Owner	delta	Points Count	65536	Pulse Sequence	qNMR_experiment.
Solvent	ACETONITRILE-d3	Spectrum Offset (Hz)	3014.5381	Sweep Width (Hz)	15007.50	Temperature (degree C)	25.000

品名
イプフルフェノキン

