

## SPECIFICATION

## Ammonium Acetate

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White crystals
Assay (CH <sub>3</sub> COONH <sub>4</sub> )	min.97.0%(mass/mass)
Solubility in water	to pass test
Residue after ignition (as sulfate)	max.0.01%(mass/mass)
pH of a 50g/l solution at 25°C	6.5~7.5
Chloride (Cl)	max.5ppm(mass/mass)
Nitrate (NO <sub>3</sub> )	max.0.001%(mass/mass)
Phosphate (PO <sub>4</sub> )	max.3ppm(mass/mass)
Sulfate (SO <sub>4</sub> )	max.0.001%(mass/mass)
Sodium (Na)	max.0.001%(mass/mass)
Potassium (K)	max.0.001%(mass/mass)
Copper (Cu)	max.2ppm(mass/mass)
Magnesium (Mg)	max.2ppm(mass/mass)
Calcium (Ca)	max.5ppm(mass/mass)
Lead (Pb)	max.3ppm(mass/mass)
Iron (Fe)	max.2ppm(mass/mass)

JIS K 8359:2006

(\*:Additional test performed by Wako but not required by JIS)