

Code No.01-0333

SPECIFICATION

Diammonium Hydrogenphosphate

Wako 1st Grade

| REQUIREMENT | SPECIFICATION |
|--------------------------------|--------------------------------------|
| Appearance | White, crystals ~ crystalline powder |
| Solubility in water | to pass test |
| pH of a 50g/l solution at 25°C | 7.7~8.3 |
| Chloride (Cl) | max.0.005% |
| Sulfate (SO ₄) | max.0.01% |
| Heavy metals (as Pb) | max.0.002% |
| Arsenic (As) | max.3ppm |
| Iron (Fe) | max.0.002% |
| Assay | min.98.5% |