

Code No.01-0449

SPECIFICATION

Antimony(III) Chloride

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay Solubility in dil. hydrochloric acid Sulfur compounds (as SO ₄) Sodium (Na) Copper (Cu) Lead (Pb) Arsenic (As) Iron (Fe)	White ~ slightly yellow, crystals ~ crystalline powder min.98.0% to pass test max.0.01% max.0.005% max.5ppm max.0.003% max.0.002% max.0.002%

JIS K 8400-1994(*:Additional test performed by Wako but not required by JIS)