

Code No.01-0481

SPECIFICATION

L-Asparagine Monohydrate
Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White, crystals ~ crystalline powder
Assay (C ₄ H ₈ N ₂ O ₃ ·H ₂ O)	min.99.0%(mass/mass)
Solubility in water	to pass test
Specific rotation [α] _{D20} (after drying at 130℃)	+33.0~+36.5°
Loss on drying at 130℃	11.5~12.5%(mass/mass)
Residue after ignition (as sulfate)	max.0.05%(mass/mass)
Chloride (Cl)	max.0.01%(mass/mass)
Sulfate (SO ₄)	max.0.02%(mass/mass)
Lead (Pb)	max.0.001%(mass/mass)
Arsenic (As)	max.1ppm(mass/mass)
Iron (Fe)	max.0.001%(mass/mass)
Ammonium (NH ₄)	max.0.02%(mass/mass)
Foreign amino acids	to pass test

JIS K 8021:2024(*:Additional test performed by Wako but not required by JIS)