

SPECIFICATION

Amino Acids Mixture Standard Solution, Type H
for Amino Acid Automated Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Cl ion concentration (Cl-)	0.10~0.13mol/l
Practical test	to pass test
L-Aspartic Acid	2.48~2.52 μ mol/ml
L-Threonine	2.48~2.52 μ mol/ml
L-Serine	2.48~2.52 μ mol/ml
L-Glutamic Acid	2.48~2.52 μ mol/ml
L(-)-Proline	2.48~2.52 μ mol/ml
Glycine	2.48~2.52 μ mol/ml
L-Alanine	2.48~2.52 μ mol/ml
L(-)-Cystine	2.48~2.52 μ mol/ml
L-Valine	2.48~2.52 μ mol/ml
L-Methionine	2.48~2.52 μ mol/ml
L-Isoleucine	2.48~2.52 μ mol/ml
L-Leucine	2.48~2.52 μ mol/ml
L-Tyrosine	2.48~2.52 μ mol/ml
L-Phenylalanine	2.48~2.52 μ mol/ml
L-Lysine	2.48~2.52 μ mol/ml
L-Histidine	2.48~2.52 μ mol/ml
Ammonium Chloride	2.48~2.52 μ mol/ml
L-Arginine	2.48~2.52 μ mol/ml