

Code No.01-0870

SPECIFICATION

Ammonium Amidosulfate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (NH ₄ OSO ₂ NH ₂) Solubility in water Residue after ignition (as sulfate) pH of a 50g/l solution at 25°C Chloride (Cl) Sulfate (SO ₄) Copper (Cu) Lead (Pb) Iron (Fe)	White crystals min.98.5%(mass/mass) to pass test max.0.01%(mass/mass) 4.7~5.7 max.0.001%(mass/mass) max.0.04%(mass/mass) max.2ppm(mass/mass) max.2ppm(mass/mass) max.5ppm(mass/mass)

JIS K 8588:2011

(*:Additional test performed by Wako but not required by JIS)