

Code No.01-1128

SPECIFICATION

Acetone for Pesticide Residue And Polychlorinated Biphenyl Tests

Concentration 300

REQUIREMENT	SPECIFICATION
* Appearance * Suitability for pesticide residue analysis (GC-FPD) Assay (CH ₃ COCH ₃) (GC) Density (20°C) Water Residue after evaporation Suitability for pesticide residue-PCB analysis	Colorless clear liquid to pass test min.99.5%(area/area) 0.789~0.792g/mL max.0.3%(mass/mass) max.2ppm(mass/mass) to pass test

JIS K 8040:2020(*:Additional test performed by Wako but not required by JIS)