

Code No.01-1627

## SPECIFICATION

### Ammonium Tetraborate Tetrahydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in water	to pass test
pH of a 10g/l solution at 25°C	8.8~9.1
Chloride (Cl)	max.0.001%
Sulfate (SO <sub>4</sub> )	max.0.005%
Heavy metals (as Pb)	max.0.001%
Iron (Fe)	max.0.001%
Assay	min.98.0%