

Code No.01-1910

SPECIFICATION

Acetonitrile

for Spectrochemical Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.780~0.782g/ml
Refractive index 20°C	1.343~1.346
Absorbance	200nm:max.0.05 210nm:max.0.03 220nm:max.0.02 225~400nm:max.0.01
Water	max. 0.01%
Fluorescence test	to pass test
Assay (capillary GC)	min.99.8%