

## SPECIFICATION

## Alisol A

for the Japanese Pharmacopoeia Crude Drugs Test (for Thin-Layer Chromatography)

REQUIREMENT	SPECIFICATION
Description (JP method)	White ~ slightly yellow, powder
Optical rotation [ $\alpha$ ]D20 (JP method)	+86~+106°
Purity (JP method)	-
Related substances (JP method)	to pass test
* Solubility in methanol (after drying)	to pass test
* Loss on drying over silica gel	max.2.0%
* Assay (HPLC)	min.95.0%

\*:Additional test performed by Wako