

SPECIFICATION

1mol/l Ammonium Dihydrogenphosphate Solution for High Performance Liquid Chromatography

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Absorbance (220nm)	max.0.10
Absorbance (230nm)	max.0.08
Absorbance (240nm)	max.0.06
Absorbance (254~400nm)	max.0.05
pH at 25°C	3.8~4.2
Fluorescence test	to pass test
Concentration (20°C)	0.95~1.05mol/l