

Code No.01-2154

SPECIFICATION

Albumin, Human, recombinant expressed in plants
for Cell Culture

| REQUIREMENT | SPECIFICATION |
|--|----------------------------------|
| Appearance | White ~ brown, crystals ~ powder |
| Solubility in water | to pass test |
| Loss on drying over silica gel in vacuum | under investigation |
| pH of a 100g/L solution at 25°C | 6.5~8.5 |
| Bacterial endotoxins | under investigation |
| Total nitrogen (N) | under investigation |
| Assay (SDS-PAGE) | under investigation |

Source:Plant expressed human serum albumin