

Code No.01-2181

SPECIFICATION

L-Aspartic Acid Standard for Amino Acid Analysis

REQUIREMENT	SPECIFICATION
Appearance	White crystalline powder
Solubility in dil.hydrochloric acid	to pass test
Specific rotation $[\alpha]_{D20}(c=8,HCl)$	$+24.9 \sim +25.8^{\circ}$
Loss on drying at 105°C	max.0.2%
Suitability for amino acids determination	to pass test
Assay (volumetric analysis)	min.99.0%
Assay (subtracting method)	measured value