

## SPECIFICATION

 $\beta$ -Alanine Standard  
for Amino Acid Analysis

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in water	to pass test
Specific rotation [ $\alpha$ ]D20(c=10,dil.HCl) (after drying)	under investigation
Loss on drying at 105°C	max.0.3%
Suitability for amino acids determination	under investigation
Assay (volumetric analysis)	min.99.0%
Assay (subtracting method)	under investigation