

## SPECIFICATION

## Amidosulfuric Acid

TraceSure(R)

REQUIREMENT	SPECIFICATION
* Appearance	White crystals
Assay (HOSO <sub>2</sub> NH <sub>2</sub> )	min.99.90%(mass/mass)
Water-insoluble matter	max.0.01%(mass/mass)
Loss on drying at 50°C	max.0.1%(mass/mass)
Residue after ignition	max.0.005%(mass/mass)
Chloride (Cl)	max.5mg/kg
Nitrate (NO <sub>3</sub> )	max.5mg/kg
Copper (Cu)	max.1 mg/kg
Lead (Pb)	max.1 mg/kg
Iron (Fe)	max.2mg/kg

JIS K 8005:2016

\*:Additional test performed by Wako but not required by JIS