

Code No.01-2377

SPECIFICATION

Arbutin

for the Japanese Pharmacopoeia Crude Drugs Test (for Assay and Thin-Layer Chroma

REQUIREMENT	SPECIFICATION
Description (JP method)	Colorless ~ white crystals or crystalline powder
Melting point (JP method)	199~201°C
Absorbance E1%1cm (280nm) (JP method)	70~76
Purity (JP method)	-
Related substances (1) (JP method)	to pass test
Related substances (2) (JP method)	to pass test
* Solubility in water	to pass test
* Loss on drying over silica gel in vacuum	max.1.0%

*:Additional test performed by Wako