

Code No.01-2401

SPECIFICATION

Anion Mixture Standard Solution 1

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Fluoride ion (F-) (20°C)	19.0~21.0mg/l
Chloride ion (Cl-) (20°C)	19.0~21.0mg/l
Nitrite ion (NO ₂ -) (20°C)	95~105mg/l
Bromide ion (Br-) (20°C)	95~105mg/l
Nitrate ion (NO ₃ -) (20°C)	95~105mg/l
Phosphate ion (PO ₄ ³⁻) (20°C)	190~210mg/l
Sulfate ion (SO ₄ ²⁻) (20°C)	95~105mg/l