

Code No.01-2814

SPECIFICATION

Amino Acids Mixture Standard Solution, Type B (High Range)

for Amino Acid Automated Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
4-Aminobutyric Acid	2.38~2.62 μ mol/mL
2-Aminoethanol	2.38~2.62 μ mol/mL
Ammonium Chloride	2.38~2.62 μ mol/mL
5-Hydroxy-DL-lysine	2.38~2.62 μ mol/mL
L-Ornithine	2.38~2.62 μ mol/mL
L-Lysine	2.38~2.62 μ mol/mL
1-Methyl-L-histidine	2.38~2.62 μ mol/mL
L-Histidine	2.38~2.62 μ mol/mL
3-Methyl-L-histidine	2.38~2.62 μ mol/mL
L-Anserine	2.38~2.62 μ mol/mL
L-Carnosine	2.38~2.62 μ mol/mL
L-Arginine	2.38~2.62 μ mol/mL