

SPECIFICATION

2-Aminoethyl Diphenylborinate

for TLC

REQUIREMENT	SPECIFICATION
Description	White ~ brown, crystals, crystalline powder or powder
Identification	to pass test
* Solubility in methanol	to pass test
* Melting point	under investigation
* Loss on drying over silica gel in vacuum	under investigation
* Assay (HPLC)	min.98.0%

based on The Japanese standards for non-Pharmacopoeial crude drugs [2-Aminoethyl Diphenylborinate](*:Additional test performed by Wako)