

Code No.01-2839

# SPECIFICATION

## Albumin, from Bovine Serum, Protease Free, Low IgG for Biochemistry

| REQUIREMENT                     | SPECIFICATION   |
|---------------------------------|---|
| Appearance                      | Slightly brown ~ brown, crystals ~ crystalline powder |
| Solubility in water             | to pass test  |
| Loss on drying at 135°C         | max.5.0%  |
| pH of a 100g/L solution at 25°C | 6.5~7.5   |
| IgG                             | under investigation                                   |
| Protease                        | to pass test  |
| Assay (electrophoresis)         | min.98%   |