

Code No.01-2847

SPECIFICATION

L-Ascorbic Acid 2-Phosphate Trisodium Salt  
Wako Special Grade

| REQUIREMENT                    | SPECIFICATION                                |
|--------------------------------|--|
| Appearance                     | White ~ pale yellowish brown, powder or mass |
| Solubility in water            | to pass test                                 |
| Water                          | under investigation                          |
| pH of a 10g/L solution at 25°C | 9.0~10.0                                     |
| Assay (HPLC)                   | min.98.0%                                    |