

## SPECIFICATION

## L-Ascorbic Acid Phosphate Magnesium Salt n-Hydrate

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly yellow, powder or mass
Identification	to pass test
Solubility in water	to pass test
Loss on drying	max.29%
Free phosphoric acid	max.1%
Heavy metals	max.20ppm
Arsenic (as As <sub>2</sub> O <sub>3</sub> )	max.4ppm
Assay (HPLC)	min.90.0%
Bacterial endotoxins	less than 10EU/g