

SPECIFICATION

Barium Acetate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay [Ba(CH ₃ COO) ₂] Solubility in water pH of a 50g/L solution at 25°C Chloride (Cl) Oxidizing substances (as NO ₃) Sodium (Na) Potassium (K) Copper (Cu) Calcium (Ca) Strontium (Sr) Lead (Pb) Iron (Fe)	White, crystals ~ crystalline powder min.99.0%(mass/mass) to pass test 7.4~8.5 max.5ppm(mass/mass) max.0.005%(mass/mass) max.0.01%(mass/mass) max.0.01%(mass/mass) max.5ppm(mass/mass) max.0.05%(mass/mass) max.0.1%(mass/mass) max.5ppm(mass/mass) max.5ppm(mass/mass)

JIS K 8376 : 2019(*:Additional test performed by Wako but not required by JIS)