

Code No.02-1276

SPECIFICATION

Bis(3-methyl-1-phenyl-5-pyrazolone)

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White ~ slightly yellow, crystalline powder ~ powder
Assay	min.97.0%
Solubility in N,N -dimethylformamide	to pass test
Residue after ignition (as sulfate)	max.0.1%
Suitability for cyanide determination	to pass test

JIS K 9545-1994

(*:Additional test performed by Wako but not required by JIS)