

Code No.02-1276

## SPECIFICATION

## Bis(3-methyl-1-phenyl-5-pyrazolone)

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay Solubility in N,N -dimethylformamide Residue after ignition (as sulfate) Suitability for cyanide determination	White ~ slightly yellow, crystalline powder ~ powder min.97.0%  to pass test  max.0.1%  to pass test

JIS K 9545-1994

(\*:Additional test performed by Wako but not required by JIS)